

Kul-PCR protokoll

Datum: 080929

Projekt: AcceLa / Nemeskodema

Primers: 28S part 1F / 1642

	µl/prov	mix x
Vatten	22	
Primer 1	0.15	
Primer 2	0.15	
BSA		
DNA	1	

Program:		File
95°C	3'	
95°C	30"	
58-53°C	30"	ic TD x1
72°C	30"	
52°C	30"	x34
72°C	30"	
72°C	10'	

Maskin nr: egenskap

Taxa:

- | | |
|---------------------|----------------------|
| 1. <u>08-85 om</u> | 11. <u>08-96 om</u> |
| 2. <u>08-86 OK</u> | 12. <u>08-98 OK</u> |
| 3. <u>08-87 om</u> | 13. <u>08-99 om</u> |
| 4. <u>08-88 om</u> | 14. <u>08-100 OK</u> |
| 5. <u>08-90 OK</u> | 15. <u>08-102 om</u> |
| 6. <u>08-91 om</u> | 16. <u>08-103 OK</u> |
| 7. <u>08-92 OK</u> | 17. <u>08-104 OK</u> |
| 8. <u>08-93 OK</u> | 18. <u>08-105 om</u> |
| 9. <u>08-94 OK</u> | 19. <u>08-106 om</u> |
| 10. <u>08-95 om</u> | 20. <u>Kontroll.</u> |

Kommentarer:

