

PCR protokoll

Datum: 081017 Projekt: Acoela/Nemertoderma.

Primers: 285 partIF / 1600R

Program:

	ul/prov	mix X
H ₂ O	22	
primer 1	0.5	
primer 2	0.5	
DNA	2	

95°C 3min	X1
94°C 30 s	X 7
58°C - 30s	12 TD
72°C 30s	
94°C 30"	x 33
52°C 30"	
72°C 30"	
72°C 10 min	X1

Maskin nr:

1	08-87	21
2	08-88	22
3	08-91	23
4	08-95	24
5	08-96	25
6	08-99	26
7	08-102	27
8	08-105	28
9	Kontroll.	29
10		30
11		31
12		32
13		33
14		34
15		35
16		36
17		37
18		38
19		39
20		40

