

Kul-PCR protokoll

Datum: 081125

Projekt: Acetab Nematoderma 2008 sept

Primers: Lco/Hco

	µl/prov	mix x
Vatten	22	
Primer 1	0.5	
Primer 2	0.5	
BSA		
DNA	2	

Program:	File
95°C 3'	
95°C 30"	
96°C 30"	
72°C 30"	
72°C 10'	

X42

Maskin nr: Experiment 2

Taxa:

- | | |
|-------------------|---------------------|
| 1. <u>08-86</u> | 11. <u>08-116</u> |
| 2. <u>08-102</u> | 12. <u>08-117</u> |
| 3. <u>08-103</u> | 13. <u>08-118</u> |
| 4. <u>08-107</u> | 14. <u>08-119</u> |
| 5. <u>08-108</u> | 15. <u>08-120</u> |
| 6. <u>08-109</u> | 16. <u>08-121</u> |
| 7. <u>08-110</u> | 17. <u>08-122</u> |
| 8. <u>08-117</u> | 18. <u>08-123</u> |
| 9. <u>08-113</u> | 19. <u>Kontroll</u> |
| 10. <u>08-115</u> | 20. |

Kommentarer:

