

# Kul-PCR protokoll

Datum: 09/005 Projekt: Acoela?

Primers: S30/SEK, 4FB/ROGR, LCO/HCO, 1200F/2450K



	µl/prov	mix x
Vatten	22	
Primer 1	0,5	
Primer 2	0,5	
BSA		
DNA	2	

Program:	File
<i>Se tidigare</i>	
<i>Program</i>	

Maskin nr: 2720 + ejmerad

Taxa:

1. 09-78 \_\_\_\_\_ 11. \_\_\_\_\_
2. \_\_\_\_\_ 12. \_\_\_\_\_
3. \_\_\_\_\_ 13. \_\_\_\_\_
4. \_\_\_\_\_ 14. \_\_\_\_\_
5. \_\_\_\_\_ 15. \_\_\_\_\_
6. \_\_\_\_\_ 16. \_\_\_\_\_
7. \_\_\_\_\_ 17. \_\_\_\_\_
8. \_\_\_\_\_ 18. \_\_\_\_\_
9. \_\_\_\_\_ 19. \_\_\_\_\_
10. \_\_\_\_\_ 20. \_\_\_\_\_



Kommentar: