

Kul-PCR protokoll

Datum: 110516

Projekt: Acacia 2010 (kort)

Primers: U178/1634L - 1200F/2450R - U2176/3449L

	µl/prov	mix x
Vatten	22	
Primer 1	0.15	
Primer 2	0.15	
BSA		
DNA	2	

Program:	File

Se tidigare
Program

Maskin nr: _____

Taxa:

- | | |
|---------------------|---------------------|
| 1. <u>U510-171</u> | 11. <u>U510-214</u> |
| 2. <u>U510-173</u> | 12. <u>U510-217</u> |
| 3. <u>U510-176</u> | 13. <u>U510-219</u> |
| 4. <u>U510-187</u> | 14. <u>Kontroll</u> |
| 5. <u>U510-192</u> | 15. _____ |
| 6. <u>U510-194</u> | 16. _____ |
| 7. <u>U510-196</u> | 17. _____ |
| 8. <u>U510-199</u> | |
| 9. <u>U510-203</u> | |
| 10. <u>U510-209</u> | |

Kommentarer:

