

PCR protokoll

Datum: 120118

Projekt: Acoela2011

Primers: TimA/TimB- S30/5FK- 4FB/1806R

	ul/prov	mix X
H ₂ O	22	
primer 1	0,5	
primer 2	0,5	
DNA	2	

Program:

95°C	4min	X1
94°C	30"	X35
50 °C	30"	
72°C	30"	
72°C	10 min	X1

Maskin nr: Geneamp2

1	11-206 Mecynostomum auritum
2	11-238 Mecynostomum sp-6
3	11-250 Mecynostomum sp-6
4	Kontroll
5	
6	
7	
8	
9	
10	
11	
12	
13	
14	
15	
16	
17	
18	
19	
20	
21	
22	
23	
24	

