

PCR protokoll

Datum: 100614-15 Projekt: Acoela

Primers: LCO/HCO - U178/1634 - UJ2176/3449L - 1200F/2450R

Program:

	ul/prov	mix X
H ₂ O	22	
primer 1	0,5	
primer 2	0,5	
DNA	2	


9 °C	min	X1
94°C	30 s	X
°C	s	
72°C	s	
72°C	min	X1

setidjau
program

Maskin nr:

	U178	COL		
1	UJ1002	2,3	2,6	21
2	UJ10-03			22
3	UJ10-10			23
4	UJ10-11			24
5	UJ10-21			25
6	UJ10-24			26
7	UJ10-29			27
8	UJ10-35			28
9	K-			29
10				30
11				31
12				32
13				33
14				34
15				35
16				36
17				37
18				38
19				39
20				40

1 2 3 4 5 6 7 8 K



1200F

UJ2176-

8 + 6 5 4 3 2 1

U178
COL