

PCR protokoll

Datum: 110817

Projekt: Acoela-2010 Karin Nilsson

Primers:

1200F/2450R - U178/1634L - UJ2176/3449L

Program:

	ul/prov	mix X	95°C 4min	95°C 4min	X1
H ₂ O	22		94°C 30" x40	94°C 30"	X2
primer 1	0,5		52°C 30"	56-54°C 30"	
primer 2	0,5		72°C 60"	72°C 60"	
DNA	2			94°C 30" 52°C 30" 72°C 60"	x36
			72°C 10 min	72°C 10 min	X1

Maskin nr:

1	UJ10-85-3Actinoposthiidae	U178/1634L
2	UJ10-196	
3	Kontroll	
4	UJ10-222	UJ2176/3449L
5	UJ10-232	
6	UJ10-235	
7	UJ10-237	
8	Kontroll	
9	UJ10-246	
10	UJ10-248	
11	UJ10-83	
12	UJ10-85-3Actinoposthiidae	
13	UJ10-97	
14	UJ10-85- sagittiferidae	
15	UJ10-243	
16	UJ10-85-3Actinoposthiidae	1200/2450R
17	Kontroll	
18		
19		
20		
21		
22		
23		
24		



Keyvan