

PCR protokoll

Datum: 120118

Projekt: Acoela 2011

Primers: U178/1634L-1200F/2450R-UJ2176/3449L

	ul/prov	mix X
H ₂ O	22	
primer 1	0,5	
primer 2	0,5	
DNA	2	

Program: 1200F/2450L

U178/1634/1642L

95°C	4min	X1
94°C	30"	X40
52°C	30"	
72°C	1'	
72°C	10 min	X1

94°C 30"
56-54°C 30" x2
72°C 1'

94°C 30"
52°C 30" x36
72°C 1'
72°C 10'

Maskin nr:

1	UJ11-01	U178/1642L
2	UJ11-17	
3	UJ11-38	
4	UJ11-47	U178/1642L
5	UJ11-67	
6	UJ11-74	
7	UJ11-84	
8	UJ11-87	
9	UJ11-89	U178/1642L
10	UJ11-98	
11	UJ11-118	
12	UJ11-134	
13	UJ11-137	
14	UJ11-145	
15	UJ11-149	
16	UJ11-169	
17	kontroll	
18		
19		
20		
21		
22		
23		
24		



U178/
1634L +
1642L

1200F/
2450L

Kayvan