

PCR protokoll

Datum: 120502

Projekt: Acoela 2011

Primers: U178/1634L - U178/28S-part1R

28S-300F/1634L


	ul/prov	mix X
H ₂ O	22	
primer 1	0,5	
primer 2	0,5	
DNA	2	

Program:

95°C	4min	X1
94°C	30"	X2
56-52°C	30"	
72°C	1'	
94C	30	34
52C	30	
72C	1min	
72°C	10 min	X1

Maskin nr:

1	UJ11-89	U178/1634L 300F/1634L
2	UJ11-118	---
3	UJ11-302	---
4	UJ11-303	---
5	UJ11-329	---
6	UJ11-379	---
7	UJ11-388	---
8	Kontroll	---
9	UJ11-159	U178/28S-part1R
10	UJ11-167	---
11	Kontroll	
12		
13	120504	
14	UJ11-225-II	U178/1634
15	UJ11-149	---
16	Kontroll	
17	UJ11-225 II	1200F/2450
18	UJ11-149	---
19	UJ11-225 II	UJ2176/3449
20	Kontroll	---
21		
22		
23		
24		



Keyvan